

CROSS PROTECTION



IN-HOUSE
MANUFACTURING
PRODUCTS
www.crossprotection.com.my



CROSS PROTECTION®



CROSS PROTECTION (M) SDN BHD
No.11, Jalan Perindustrian 5, Batu 5 1/2,
Off Jalan Haji Abdul Manan, 41050 Klang,
Selangor, Malaysia.



operation@crossprotection.com /
admin@crossprotection.com



+603 33920000 / +6016 2620436

GLOVES

Medical gloves are **Disposable Gloves** used during medical examinations and procedures to help prevent cross-contamination between caregivers and patients.



XPERT LATEX POWDER EXAM GLOVE

The superior protection will keep safe from exposing to high risk task in longer duration.

Available sizes : **XS S M L XL**



XTEND LATEX POWDER EXAM GLOVE

Mid-level protection, balanced with excellent flexibility and dexterity for lighter-duty applications.

Available sizes : **XS S M L XL**



XCEL LATEX POWDER FREE EXAM GLOVE

Delivering an amazing fit and level of comfort. It also specialized formulated tacky grip.

Available sizes : **XS S M L XL**



XCED LATEX POWDER FREE EXAM GLOVE

Value general-purpose latex disposable gloves that provides fair barrier protection and chemical resistance

Available sizes : **XS S M L XL**



XTRETCH NITRILE EXAM GLOVES

Hand protection with superior strength. Great comfort and tactile sensitivity.

Available sizes : **XS S M L XL**



COLOR LATEX POWDER FREE GLOVES

Design with perfect smell while keeping you protected from chemicals and bacterial.

Available sizes : **XS S M L XL**

Available colors :



XTRETCH BASIC POWDER FREE EXAM GLOVES

Lightweight version of popular Xtretch nitrile gloves. Ambidextrous and suitable for general usage.

Available sizes : **XS S M L XL**

Available colors :



XTREME DIAMOND POWDER FREE GLOVES

Embossed diamond texture on front and back with exceptional grip

Available sizes : **M L XL XXL**

Available colors :



FACE MASK

A surgical **Face Mask** is a disposable medical device that protects from bacteria shed in liquid droplets and aerosols from the wearer's mouth and nose.



RES-Q 300 3 Ply Surgical Face Mask

It provides you the dependable protection against fluid splashes at **80mmHg** while remaining excellent breathability.

Available in : **Earloop & Tie-On**

Protection level : **ASTM LEVEL 1 / EN TYPE II**

Available color :



RES-Q 300 PLUS 3 Ply Surgical Face Mask

The highly efficient filtration media. Resist moderate fluid splashes at **120mmHg** while maintaining smooth breathability.

Available in : **Earloop & Tie-On**

Protection level : **ASTM LEVEL 2 / EN TYPE IIR**

Available color :



RES-Q 300 ULTRA 3 Ply Surgical Face Mask

Excellent surgical mask for those who need highest level of protection. Tested and passed splash resistant **160mmHg**.

Available in : **Earloop & Tie-On**

Protection level : **ASTM LEVEL 3 / EN TYPE IIR**

Available color :



ANTIFOG SHIELD 3 Ply Surgical Face Mask

Extra protection with excellent surgical mask for medical practitioners. Tested and passed splash resistant **120mmHg**.

Available in : **Earloop & Tie-On**

Protection level : **ASTM LEVEL 2 / EN TYPE IIR**

Available color :



CHILDREN FACE MASK 3 Ply Surgical Face Mask

Our Children Face Mask is designed especially for children from age of 4 to 12. preventing child from inhaling small particles and bacteria.

Available in : **Earloop**

Protection level : **ASTM LEVEL 2 / EN TYPE IIR**

Available color :



DENTAL ACCESSORIES



DENTAL BIB

- Made of 2 ply paper and 1 ply PE
- Extra heavy weight 18+18gsm origin paper
- 45cm x 33cm size paper
- Excellent absorbency
- Hygiene & Comfortable Personal Barrier
- Available color



SALIVA EJECTOR 145mm

- Dimension 145mm
- PVC, coated steel wire
- Cylindrical form
- Bend easily
- Smooth edge and suction, no block
- Assorted color

PROTECTION APPAREL

AAMI SURGEON CAP

Material : Premium SMMMMS
Standard : Easy tie, air breathability,
Water repellent, anti-static
Color : 



ROUND CAP

Material : SMS
Type : Basic / Premium
Standard : Stretchable
Color : 



ROUND CAP (LATEX FREE)

Material : Premium SMMMMS
Standard : Antistatic performance
Spray and water resistant
Color : 



SHOE COVER

Material : PP spunbond
Type : Standard / Antiskid
Color : 



SLEEVE COVER

Material : Premium SMMMMS
Color : 

ISOLATION GOWN

Material : Premium SMMMS
Standard : Splash resistant, waist ties
Color : 



SURGICAL GOWN

Standard : AAMI level 1, 2, 3
Splash and fluid protection
Soft breathable
Non- irritating material
Low-linting
Type : Standard / Reinforced
Content : 1x Gown
2x Paper towel
1x Outer wrap
Size : S M L XL
Color : 

STERILIZATION PRODUCTS

Packing the medical devices before sterilization and maintain the sterility of the products after sterilization.
With triple side seal for secure closure and double indicator (inside and outside of the pouch). EO and Steam duo indicator. CE certified.



POUCH

Duo indicator



SIZES	PACKING
60mm X 130mm	200 pcs X 30 boxes
70mm X 260mm	200 pcs X 40 boxes
90mm X 260mm	200 pcs X 20 boxes
135mm X 280mm	200 pcs X 20 boxes
190mm X 360mm	200 pcs X 8 boxes
300mm X 490mm	200 pcs X 5 boxes



FLAT REEL

Tri indicator



SIZES	PACKING
50mm X 200m	16 rolls per carton
75mm X 200m	10 rolls per carton
100mm X 200m	6 rolls per carton
150mm X 200m	4 rolls per carton
200mm X 200m	4 rolls per carton
250mm X 200m	2 rolls per carton
300mm X 200m	2 rolls per carton
350mm X 200m	2 rolls per carton
400mm X 200m	2 rolls per carton



GUSSETED REEL

Tri indicator



SIZES	PACKING
75mm X 38mm X 100m	8 rolls per carton
100mm X 50mm X 100m	6 rolls per carton
150mm X 55mm X 100m	4 rolls per carton
200mm X 65mm X 100m	2 rolls per carton
250mm X 70mm X 100m	2 rolls per carton
300mm X 80mm X 100m	2 rolls per carton
400mm X 100mm X 100mm	2 rolls per carton



TYVEK REEL

SIZES	PACKING
75mm X 70m	8 rolls per carton
100mm X 70m	6 rolls per carton
150mm X 70m	4 rolls per carton
200mm X 70m	4 rolls per carton
250mm X 70m	2 rolls per carton
300mm X 70m	2 rolls per carton



STERILIZATION WRAPPER

Material : SMS

SIZES	PACKING
50cm X 50cm	500pcs per carton
60cm X 60cm	500pcs per carton
75cm X 75cm	250pcs per carton
90cm X 90cm	250pcs per carton
100cm X 100cm	250pcs per carton
120cm X 120cm	125pcs per carton
120cm X 140cm	125pcs per carton
120cm X 160cm	125pcs per carton
120cm X 180cm	125pcs per carton



GREEN CREEP PAPER

SIZES	PACKING
60cm X 60cm	500pcs per carton
75cm X 75cm	250pcs per carton
90cm X 90cm	250pcs per carton
100cm X 100cm	250pcs per carton
120cm X 120cm	125pcs per carton
120cm X 140cm	125pcs per carton
120cm X 160cm	125pcs per carton
120cm X 180cm	125pcs per carton