



CROSS PROTECTION®



2019

CATALOG

Manufacturer of Dental and
Medical Disposable Products



CROSS PROTECTION®



CROSS PROTECTION, established in 1997, 100% Malaysian owned manufacturing company, located in Klang, Selangor. We are having 5 main production divisions : Sewing, Plastic Injection and Extrusion, Latex and Nitrile Gloves, Dental Bib, and Surgical Mask.

Our existing inhouse production included cleanroom surgical mask, surgical gown, cap, shoe cover, patient bedsheet and pillow case, gown, saliva ejector, dental bib, latex and nitrile gloves, surgical sets, sterilization drapes and customized procedure sets like caesarean, ophthalmic and etc. All products are registered in Malaysia and Europe under MRHA (UK) and MCCAA (EU).

With continuous research and developments, every piece of products is push for the perfection finish. All products are manufactured in clean room environments. In CROSS PROTECTION, we believe **QUALITY** should come first in Medical Industry. Every medical practitioner deserve for the best protection when they are helping the patient.

We strive to provide top protection and quality products to you and to become your valuable partners during your professional path.

Gloves

ASTM RATED DISPOSABLE GLOVES

- PAGE** 5 • XPERT Latex Powdered - Advance Stretch Extra Durable
 • XTEND Latex Powdered - Standard
 • XCEL Latex Powder Free - Tacky Grip Extra Durable
 • XCED Latex Powder Free - Standard
 • COLOR Latex Powder Free - Color Gloves
- 6 • XTRETCH BASIC Blue Nitrile Glove - Powder Free
 • XTRETCH Violet Blue Nitrile Glove - Powder Free
 • XTREME Powder Free - Multipurpose Nitrile Gloves
 • XTREME (pack) Powder Free - Multipurpose Nitrile Gloves

Gloves

NON ASTM RATED DISPOSABLE GLOVES

- PAGE** 8 • Clear Color Plastic Gloves
 • Clear Color Plastic Shoulder Gloves
 • Yellow Color Utility Gloves
 • Yellow Color Oven Gloves
 • Silver Color Oven Gloves
 • Yellow Color Cold Resistant Gloves

Face Mask

ASTM RATED DISPOSABLE FACE MASK

- PAGE** 10 • Res-Q 300 (ASTM Level 1) - 3ply Surgical Face Mask
 • Res-Q 300 Plus - (ASTM Level 2) - 3ply Surgical Face Mask
 • Res-Q 300 Ultra - (ASTM Level 3) - 3ply Surgical Face Mask
 • Antifog Shield - (ASTM Level 2) - 3ply Surgical Face Mask
- 11 • 3ply Surgical Face Mask (ASTM Level 1)
 • YouSafe Adult (ASTM Level 2) - 3ply Surgical Face Mask
 • YouSafe Children (ASTM Level 1) - 3ply Surgical Face Mask
 • YouSafe Children (ASTM Level 2) - 3ply Surgical Face Mask

Face Mask

NON ASTM RATED DISPOSABLE FACE MASK

- PAGE** 12 • Paper Face Mask
 • Face Mask - Standard

Apparel

ASTM RATED PROTECTION APPAREL

- PAGE** 15 • Sontara Gown
 • Surgical Gown - AAMI Level 2
 • Surgical Gown - AAMI Level 3
- 16 • AAMI Round Cap
 • AAMI Surgeon Cap
 • AAMI Shoe Cover
 • AAMI Isolation Gown - Level 2
 • Coverall Suit with Hood - Non Sterile

Apparel

NON ASTM RATED PROTECTION APPAREL

- PAGE** 17 • Clip Cap
 • Round Cap
 • Beard Cover
 • Astronaut Cap
- 18 • Roll Apron
 • Laboratory Coat
 • Heavyweight Cystoscopy Gown with Thumb Loop
 • Plastic Apron HDPE
 • Isolation Gown - Non Sterile
- 19 • Sleeve Cover - Non Woven with Knitted Cuff
 • Plastic Sleeve Cover - Elastic Cuff
 • Shoe Cover
 • Shoe Cover Antiskid Base
 • Boot Cover
- 20 • Plastic Shoe Cover
 • Shoe Cover Refill - for Mac Shoe Machine
 • Shoe Cover Machine

Sterilization

- PAGE** 22 • Pouches
 • Gussetted Reel
 • Flat Reel
 • Tyvek Reel
- 23 • Autoclave Tape
 • Green Crepe Paper
 • Wrap

Bandage & Tape

- PAGE** 24 • Crepe Bandage
 • Surgical Tape Paper
 • Bandage Scissors

Gauze & Sponge

- PAGE** 25 • Non Woven Sponge Roll 100m
 • Dental Roll
 • Cotton Wool Roll
 • Gauze Roll 100yrd
- 26 • Cotton Wool Ball 150gsm
 • Combined Dressing Roll
 • Precut Non Woven Sponge 4ply - Non sterile
 • Precut Gauze Sponge 8ply - Non Sterile

Hospital Admission

- PAGE** 29 • Mini Compressed Towel - White
 • Compressed Towel - Assorted Color
 • Basic Dressing Set
 • AAMI Disposable Pillow Case
 • AAMI Disposable Bedsheet
- 30 • Cardiology Stethoscope
 • SAP Absorbent Underpad
 • Cotton soft - Pillow case
 • Cotton soft - Bedsheet
- 31 • Procedure set - Cardiovascular Drape
 • Procedure set - Caesarean
 • Procedure set - Ophthalmic

Dental Accessories

- PAGE** 32 • Dental Mirror
 • Dental Probe
 • Microtip Surgical Aspirator - Sterile
 • Denture Box
 • Dappen Dish
- 34 • Saliva Ejector
 • Dental Tubing Sleeve
 • Dental Bib
 • Barrier Film Dispenser

Consumable Products

- PAGE** 35 • Paper Hand Towel
 • Hand Roll Towel
 • Protective Goggles - Clear Antifog
- 36 • Paper Hand Towel Dispenser
 • Hand Roll Towel Dispenser - Autocut
 • Waste Bag
- 37 • First Aid Bag - Small Compact
 • First Aid Bag - Large

TOP
PICKS!



CROSS PROTECTION®

COLORTM LATEX GLOVES Powder Free



with fragrance



**TOP
PICKS!**



CROSS PROTECTION®

Latex Examination Gloves





Advance Stretch Extra Durable - Powdered

This heavy-duty disposable glove will keep you safe from exposing to high risk task in longer duration.

Code	550-CP	650-CP	750-CP	850-CP	950-CP
Size	X-Small	Small	Medium	Large	X-Large

Packing: 100pcs per box, 10 boxes per carton

Material: Natural Rubber Latex

Features: • Smooth Surface • Beaded Cuff • Lightly Powdered • Ambidextrous

Color: ☐



Standard - Powdered

Mid-level protection, balanced with excellent flexibility and dexterity for lighter-duty applications.

Code	550-ECOP	650-ECOP	750-ECOP	850-ECOP	950-ECOP
Size	X-Small	Small	Medium	Large	X-Large

Packing: 100pcs per box, 10 boxes per carton

Material: Natural Rubber Latex

Features: • Smooth Surface • Beaded Cuff • Lightly Powdered • Ambidextrous

Color: ☐



Tacky Grip Extra Durable - Powder Free

It specialized formulated tacky grip will make sure nothing slips away from your hand.

Code	550-CPF	650-CPF	750-CPF	850-CPF	950-CPF
Size	X-Small	Small	Medium	Large	X-Large

Packing: 100pcs per box, 10 boxes per carton

Material: Natural Rubber Latex

Features: • Full Palm Tacky Textured • Advanced Stretch Formulation • Beaded Cuff • Ambidextrous

Color: ☐



Standard - Powder Free

Value leader general-purpose latex disposable gloves that provides fair barrier protection

Code	550-ECOCPP	650-ECOCPP	750-ECOCPP	850-ECOCPP	950-ECOCPP
Size	X-Small	Small	Medium	Large	X-Large

Packing: 100pcs per box, 10 boxes per carton

Material: Natural Rubber Latex

Features: • Micro-Roughened Surface • Advanced Stretch Formulation • Beaded Cuff • Ambidextrous

Color: ☐



with Fragrance - Powder Free

Designed with perfect while keeping you protected from chemicals and bacterial.

Code	550-C	650-C	750-C	850-C	950-C
Size	X-Small	Small	Medium	Large	X-Large

Packing: 100pcs per box, 10 boxes per carton

Material: Natural Rubber Latex

Features: • Beaded Cuff • Ambidextrous • Advanced Stretch Formulation • Full Palm Tacky Textured

Color:



Xtretch™ Basic



Blue Nitrile Gloves - Powder Free

Lightweight version of popular Xtretch nitrile gloves. Ambidextrous and suitable for general lab, clinic examination, non-hazardous usage.

Code	550-NF2	650-NF2	750-NF2	850-NF2	950-NF2
Size	X-Small	Small	Medium	Large	X-Large

Packing: 250pcs per box, 4 boxes per carton **Material:** Acrylonitrile-butadiene

Features: • Finger-Tip Textured • Beaded Cuff • Chlorinated
• Ambidextrous • Latex Free

Color:

Xtretch™



Violet Blue Nitrile Gloves - Powder Free

It offers 3 times the puncture resistance of latex, with great comfort and tactile sensitivity.

Code	550-CPNF	650-CPNF	750-CPNF	850-CPNF	950-CPNF
Size	X-Small	Small	Medium	Large	X-Large

Packing: 100pcs per box, 10 boxes per carton **Material:** Acrylonitrile-butadiene

Features: • Finger-Tip Textured • Beaded Cuff • Chlorinated
• Ambidextrous • Latex Free

Color:

Xtreme™



Multipurpose Nitrile Gloves - Powder Free

Double-sided diamond embossing provides excellent grip in wet and dry while retaining exceptional tactile sensitivity.

BLACK	Code	750-MN	850-MN	950-MN	000-MN
	Size	Medium	Large	X-Large	2X-Large

ORANGE	Code	750-XN	850-XN	950-XN	000-XN
	Size	Medium	Large	X-Large	2X-Large

Packing: 50pcs per box, 10 boxes per carton **Material:** Nitrile Rubber

Features: • Exceptional Grip • Embossed Diamond Texture • Ambidextrous
• Latex Free

Color:

Our Xtreme Multipurpose Nitrile Gloves now comes with pack. Excellent dexterity and elasticity while still being very strong and durable.

BLACK	Code	750-MN2	850-MN2	950-MN2	000-MN2
	Size	Medium	Large	X-Large	2X-Large

ORANGE	Code	750-XN2	850-XN2	950-XN2	000-XN2
	Size	Medium	Large	X-Large	2X-Large

Packing: 2pcs per pack, 50 packs per carton **Material:** Nitrile Rubber

Features: • Exceptional Grip • Embossed Diamond Texture • Ambidextrous
• Latex Free

Color:

TOP PICKS!



CROSS PROTECTION®

Xtretch™

Violet Blue Nitrile Gloves
Powder Free



TOP PICKS!



CROSS PROTECTION®

Xtreme™

**3D DIAMOND
TEXTURED NITRILE**

POWDER FREE MULTIPURPOSE GLOVES



**with
Xtreme
Grip**



Clear Color Plastic Gloves



Use these gloves for food handling, cleaning and in any situation where you do not want to dirty your hands.

Code	L0012-S	L0012-M	L0012-L
Size	Small	Medium	Large

Color: Clear

Packing: 100pcs per box, 100 bags per carton

Material: Low Density Polyethylene

Features: • Easy to slip on and off • Smooth for improved feel
• 100% latex-free • Fit either hand • Disposable

Clear Color Plastic Shoulder Gloves



Special design up to shoulder for a higher level of protection.

Code	L83003-M	L83003-L
Size	Medium	Large

Color: Clear

Packing: 100pcs per box, 10 bags per carton

Material: Low Density Polyethylene

Features: • Easy to slip on and off • Smooth for improved feel
• 100% latex-free • Fit either hand • Disposable

Yellow Color Utility Gloves



These UTILITY gloves also known as HOUSEHOLD gloves. It protects your hands during cleaning or dishwashing.

Code	650-UG	750-UG	850-UG	950-UG
Size	Small	Medium	Large	X-Large

Packing: 1 pair per bag, 144 bags per carton

Features: • Heat resistant

Color:

Yellow Color Oven Gloves



These OVEN gloves are the perfect choice when removing hot dishes from the oven and hob that can withstand heat of up to 400°C.

Code	750-HR
Size	Medium

Packing: 1 pair per bag

Color:

Features: • Heat resistant • Fit for both right and left hand

Silver Color Oven Gloves



These OVEN gloves are special designed for using in the kitchen and other high temperature work place that withstand heat of up to 400°C.

Code	775-HR
Size	Medium

Packing: 1 pair per bag

Color:

Features: • Higher heat resistant • Fit for both right and left hand

Yellow Color Cold Resistant Gloves

These COLD RESISTANT gloves have a good thermal insulation and cold resistance down to -40°C.

Code	950-FS
Size	Large

Packing: 1 pair per bag

Color:

Features: • High resistance to cold, cut, abrasion and puncture



TOP PICKS!



CROSS PROTECTION®



Res-Q 300 | ASTM Level 1

3ply Surgical Face Mask

Proven to control the spread of germs or keeping airborne dust from becoming a nuisance.

Packing: 50pcs per box, 20 boxes per carton

Features: • ASTM LEVEL 1

- Low Fluid Resistance: 80mmHg
- Bacterial Filtration Efficiency: >98%
- Particle Filtration Efficiency: >98%
- Breathability - Delta P: <3.0mmH2O/cm²
- Flammability: Class 1

Color:

EAR LOOP

Code	2-400EL
------	---------

TIE ON

Code	2-400TO
------	---------

Res-Q 300 Plus | ASTM Level 2

3ply Surgical Face Mask

The highly efficient filtration media can filter or trap all the common bacteria and particle surrounding us.

Packing: 50pcs per box, 20 boxes per carton

Features: • ASTM LEVEL 2

- Moderate Fluid Resistance: 120mmHg
- Bacterial Filtration Efficiency: >99%
- Particle Filtration Efficiency: >99% @ 0.1 micron
- Breathability - Delta P: <3.0mmH2O/cm²
- Flammability: Class 1

Color:

EAR LOOP

Code	2-406EL
------	---------

TIE ON

Code	2-406TO
------	---------

Res-Q 300 Ultra | ASTM Level 3

3ply Surgical Face Mask

Most premium Surgical mask, excellent surgical mask for those who need highest level of protection.

Packing: 40pcs per box, 20 boxes per carton

Features: • ASTM LEVEL 3

- High Fluid Resistance: 160mmHg
- Bacterial Filtration Efficiency: >99%
- Particle Filtration Efficiency: >99% @ 0.1 micron
- Breathability - Delta P: <3.0mmH2O/cm²
- Flammability: Class 1

Color:

EAR LOOP

Code	2-407EL
------	---------

TIE ON

Code	2-407TO
------	---------

Antifog Shield | ASTM Level 2

3ply Surgical Face Mask

Antifog Clear Shield for extra protection. Excellent surgical mask for medical practitioners.

Packing: 25pcs per box, 10 boxes per carton

Features: • ASTM LEVEL 2

- Moderate Fluid Resistance: 120mmHg
- Bacterial Filtration Efficiency: >99%
- Particle Filtration Efficiency: >99% @ 0.1 micron
- Breathability - Delta P: <3.0mmH2O/cm²
- Flammability: Class 1

Color:

EAR LOOP

Code	2-401EL
------	---------

TIE ON

Code	2-401TO
------	---------

Conforms and tested by USA Nelson Lab. Manufactured in ISO 7 clean room environment.

Application: Ideal for procedures where low amounts of fluid, spray or aerosols are produced.

Material: Inner (PP Spunbond), Middle (PP Meltblown Non Woven Fabric), Outer (PP Spunbond Non Woven Fabric)

3Ply Surgical Face Mask

ASTM Level 1

It provides you the dependable protection against low amounts of fluid, spray or aerosols are produced.

Packing: 3pcs per pack, 50 packs per carton

Features: • ASTM LEVEL 1

- Low Fluid Resistance: 80mmHg
- Bacterial Filtration Efficiency: >98%
- Particle Filtration Efficiency: >98%
- Breathability - Delta P: <3.0mmH2O/cm²
- Flammability: Class I

Color:

EAR LOOP

Code	2-433EL
------	---------

TIE ON

Code	2-433TO
------	---------

YouSafe - Adult | ASTM Level 2

3ply Surgical Face Mask

Adult Surgical Face Mask is designed to resist moderate fluid splashes while maintaining smooth breathability. The highly efficient filtration media can filter all the common bacteria and particle surrounding us.

Packing: 5pcs per pack, 50 packs per carton

Features: • ASTM LEVEL 2

- Moderate Fluid Resistance: 120mmHg
- Bacterial Filtration Efficiency: >99%
- Particle Filtration Efficiency: >99% @ 0.1 micron
- Breathability - Delta P: <3.0mmH2O/cm²
- Flammability: Class I

Color:

EAR LOOP

Code	2-578EL
------	---------

YouSafe - Children

3ply Surgical Face Mask

Children Face Mask is designed especially for children from age of 4 to 12. The light weight design ensures the wearer can breathe easily. Manufactured in ISO 7 clean room environment.

ASTM Level 1

Packing: 50pcs per box, 20 boxes per carton

Features: • ASTM LEVEL 1

- Low Fluid Resistance: 80mmHg
- Bacterial Filtration Efficiency: >98%
- Particle Filtration Efficiency: >98%
- Breathability - Delta P: <3.0mmH2O/cm²
- Flammability: Class I

Color:

EAR LOOP

Code	2-68BEL
------	---------

ASTM Level 2

Packing: 3pcs per pack, 50 packs per carton

Features: • ASTM LEVEL 2

- Moderate Fluid Resistance: 120mmHg
- Bacterial Filtration Efficiency: >99%
- Particle Filtration Efficiency: >99% @ 0.1 micron
- Breathability - Delta P: <3.0mmH2O/cm²
- Flammability: Class I

Color:

EAR LOOP

Code	2-678EL
------	---------

Conforms and tested by USA Nelson Lab. Manufactured in ISO 7 clean room environment.

Application: Ideal for procedures where low amounts of fluid, spray or aerosols are produced.

Material: Inner (PP Spunbond), Middle (PP Meltblown Non Woven Fabric), Outer (PP Spunbond Non Woven Fabric)

Paper Face Mask

This Paper Face Mask is made of origin paper pulp. The extra soft cotton loops is designed for maximum comfort. It is ideal for general usage but not for medical.

1ply **Code** 5-900

2ply **Code** 5-901

Packing: 100pcs per pack, 100 packs per carton

Type: Earloop

Color: 

Face Mask - Standard



3ply **Code** 2-40-EL

Filtration level up to 3 micron and it is comfortable and breathable during wearing. These face mask are also best used in food service and cleaning.

Packing: 50pcs per pack, 40 packs per carton

Type: Earloop

Features: • Fiberglass Free • 3.0 Micron Filtration • Plastic Coated Wire Nose Piece

Color:      



FACE MASK
NON-ASTM RATED
DISPOSABLE GLOVES



ARE YOU PROTECTED?

Choosing the appropriate face mask is a key component in minimizing the spread of potentially infectious diseases

**TOP
PICKS!**



CROSS PROTECTION®

SURGICAL GOWN

TOP
PICKS!



CROSS PROTECTION®

Sontara gown

**PREMIUM
SELECTION**



**Maximum
Comfort**



Air Permeability



Low Lint



High Opacity



Protection



Soft

Sontara® fabrics provide medical apparel that is highly comfortable to wear and fulfills international product norms, minimizing the risk of infections.

The remarkable softness and high drapeability of Sontara® fabrics help optimize surgeons comfort and dexterity throughout operations. Sontara® materials are opaque and minimize distracting glare. Sontara® drapes are strong and tear-resistant, yet conform easily and stay in place without slipping or sliding.

Material fully import from USA and tailored made to your size.

Sontara Gown

Where Performance and Efficiency is your main concern, our Sontara Gown is your perfect choice. **Sontara®** has a long history of manufacturing products for highly regulated end users (e.g. cleanroom, aerospace, health care and pharmaceutical).

Sontara® quality is consistent and uncompromising. The excellent quality material assures you a perfect protection during surgery. Material fully import from USA and tailored made to your size. Made in Malaysia.

Code	ST0001	ST0002	ST0003	ST0004
Size	Small	Medium	Large	X-Large

Packing: 1pc per pack, 25 packs per carton

Content: • 1pc of gown • 2pcs of paper towel • 1 outer wrap

Color: 

Surgical Gown | AAMI Level 2

Medium Weight Multi-ply Blue and Yellow Surgical Gowns

BASIC PLUS

Meets ANSI / AAMI PB70 Level 2, where risk of exposure to fluid, splash and pressure on gown are LOW.

Code	2-813	2-823	2-833	2-843
Size	Small	Medium	Large	X-Large

Packing: 1pc per pack, 25 packs per carton

Content: • 1pc of gown • 2pcs of paper towel • 1 outer wrap

Color:  

BASIC PLUS - Reinforced

With additional PE protection on both sleeves and front chest to provide extra barrier protection to potentially infectious blood and other bodily fluids.

Code	2-913	2-923	2-933	2-943
Size	Small	Medium	Large	X-Large

Packing: 1pc per pack, 25 packs per carton

Content: • 1pc of gown • 2pcs of paper towel • 1 outer wrap

Color:  

ANSI / AAMI PB70	AAMI LEVEL 2
TEST METHOD	Water resistance: Impact penetration AATCC 42 Water resistance hydrostatic pressure AATCC 127
REQUIREMENT	Spray Impact $\leq 1.0g$ Hydrostatic Pressure $\geq 20cm$

Surgical Gown | AAMI Level 3

Heavy Weight Multi-ply Blue and Yellow Surgical Gowns

PREMIUM

Meets ANSI / AAMI PB70 Level 3, where risk of exposure to fluid, splash and pressure on gown are MODERATE.

Code	2-003	2-004	2-005	2-006
Size	Small	Medium	Large	X-Large

Packing: 1pc per pack, 25 packs per carton

Content: • 1pc of gown • 2pcs of paper towel • 1 outer wrap

Color:  

PREMIUM - Reinforced

With additional PE protection on both sleeves and front chest to provide extra barrier protection to potentially infectious blood and other bodily fluids.

Code	2-803	2-804	2-805	2-806
Size	Small	Medium	Large	X-Large

Packing: 1pc per pack, 25 packs per carton

Content: • 1pc of gown • 2pcs of paper towel • 1 outer wrap

Color:  


ANSI / AAMI PB70	AAMI LEVEL 3
TEST METHOD	Water resistance: Impact penetration AATCC 42 Water resistance hydrostatic pressure AATCC 127
REQUIREMENT	Spray Impact $\leq 1.0g$ Hydrostatic Pressure $\geq 50cm$

AATCC 42 : Measure the resistance of fabrics to the liquid penetration of water by Impact
AATCC 127 : Measure the resistance of fabrics to the liquid penetration of water by Impact under constant and increasing hydrostatic pressure



ASTM RATED
DISPOSABLE GLOVES
APPAREL

	SIZING	Tolerance +/- 10%	
SIZE	GOWN LENGTH	SLEEVE LENGTH	GOWN WIDTH
S	110cm	55cm	142cm
M	120cm	57cm	150cm
L	130cm	60cm	158cm
XL	140cm	65cm	165cm



AAMI Round Cap

AAMI Round Cap with Antistatic Performance, Spray and Water Resistant, AAMI Level 2, Latex Free.

Code 2-704

Packing: 50pcs per pack, 20 packs per carton

Color: 

AAMI Surgeon Cap

Our SMMMMS material offer cotton soft feeling and greatest comfort to wearer. Every piece of surgeon cap are tailored made to suit your profession.

Code 2-703

Packing: 100pcs per pack, 10 packs per carton

Type: Tie  Elastic 

Features: • Water repellent • Air breathability • Anti-static

AAMI Shoe Cover

AAMI Shoe Cover with Antistatic Performance, Spray and Water Resistant, AAMI Level 2, Latex Free.

Code 2-888

Packing: 100pcs per pack, 10 packs per carton

Color: 

AAMI Isolation Gown | Level 2

Medium Weight Multi-ply Blue and Yellow Isolation Gowns

Comfortable, strong and durable SMS material for demanding situations.

Code 2-908

Color: 

Packing: 10pcs per pack, 10 packs per carton

Coverall Suit with Hood

NON - STERILE

STANDARD

Our Coverall Suit provides full upper body coverage and is made of non woven material. These zipper front coats are lightweight and comfortable to wear.

Code 2-902

Color: 

Packing: 5pcs per pack, 10 packs per carton


PE COATED

Our PE Coated Coverall Suit provides full upper body coverage and is made of non woven material. These zipper front coats are lightweight and comfortable to wear.

Code 2-9028

Color: 

Packing: 5pcs per pack, 10 packs per carton



Clip Cap

BASIC

This Clip Cap is made of non-woven material, breathable, lightweight, non-absorbent and economical. It is compactly folded in strip and easily open to shape.

Code 2302102

Packing: 100pcs per pack, 20 packs per carton

Color: ☐

PREMIUM

This Clip Cap is made of high quality non-woven material, breathable, lightweight and non-absorbent. It is compactly folded in strip and easily open to shape.

Code 2-702

Packing: 100pcs per pack, 10 packs per carton

Color: ☐ ☒



Round Cap

BASIC

Our Basic Round Cap is made of PP spunbond non woven fabric material. It is stretchable and comfortable during wearing.

Code 2-700

Packing: 100pcs per pack, 10 packs per carton

Color: ☐ ☒ ☒

PREMIUM

Our Round Cap Premium is made of PP spunbond non woven SMS material. It is stretchable and comfortable during wearing.

Code 2-701

Packing: 100pcs per pack, 10 packs per carton

Color: ☐ ☒ ☒



Beard Cover

This Beard Cover is ideal for use in food plants, pharmaceutical and lab environments. Full coverage for reliable contamination control.

Code 2-767

Packing: 100pcs per pack, 10 packs per carton

Color: ☐



Astronaut Cap

This Astronaut Cap is commonly found in food industry, laboratories, pharmaceutical, hospital and other various fields. It is breathable, lightweight and skid proof.

Code 2-707

Packing: 100pcs per pack, 10 packs per carton

Color: ☒



Roll Apron

Made from polythene, these aprons are packaged in a roll and dispense simply and continuously. Disposable-Anti-bacterial-Ties to secure at the back and looped around the neck.

Add on the Stainless Steel Dispenser for best organizing your space.

Code 4-208

Size: 28x55inch

Color: ☐

Thickness: 0.035

Surface: Embossed

Packing: 200pcs per roll, 5 rolls per carton

Code 4-208D

Product: Stainless Steel Dispenser

Color: ☐

Packing: Per unit



Laboratory Coat

Single Use With 3 Pockets

This Disposable Laboratory Coat is made of 35g light weight and high breathability polypropylene material. Mid length up to knee with elastic cuff, collar, pocket and snap front closures. It is autoclavable and ideal for laboratory use.

Code 4-201C

Packing: 5pcs per pack, 10 packs per carton

Color: ☐

Heavyweight Cystoscopy Gown

With Thumb Loop

"Cystoscopy Gown" provide protection for areas where heavy fluids are present. Made of extra thick polyethylene material. Unique design, thumb loop on sleeve to make donning gloves easier.

Code 4-600

Packing: 1pc per pack, 100 packs per carton

Color: ☒



Plastic Apron HDPE

This High Density Disposable Plastic Apron is made of polyethylene material. Full length cutting offers added protection for hygienic and safety in various environments.

Code 4-211

Packing: 100pcs per bag, 10 bags per carton

Type: Smooth surface

Color: ☐

Isolation Gown - Non Sterile

Made of high breathability SMS material for best comfort wearing. Knitted cuff for secure cuff closure. Registered under UK MHRA - Personal Protective Equipment.

STANDARD

Code 2-9008

Packing: 10pcs per bag, 10 bags per carton

Color: ☒

PE COATED

Code 2-9018

Packing: 10pcs per bag, 10 bags per carton

Color: ☒



Sleeve Cover

Non Woven With Knitted Cuff 1 End

This Sleeve Cover is made of non woven material. It can prevent and isolate dust, particle, blood, bacteria and virus from invading.

Code 2-708

Packing: 20pcs per pack, 10 packs per carton

Color: ●



Plastic Sleeve Cover

Elastic Cuff

Our Plastic Sleeve Cover is single use sleeve protectors, lightweight for comfort and dexterity. It is waterproof for protection against spills, splashes, dirt or residue from daily jobs.

Code 4-500

Packing: 100pcs per pack, 10 packs per carton

Color: ●



Shoe Cover

This Shoe Cover made of PP spunbond material. May be used to protect against dirt, grime and certain dry particulates in non-hazardous environments.

Code 2-865

Packing: 100pcs per pack, 10 packs per carton

Color: ●



Shoe Cover Antiskid Base

This Antiskid Base Shoe Cover offers improved floor traction. It is ideal for industrial, food/processing, construction and certain environmental cleanup applications.

Code 2-800

Packing: 100pcs per pack, 10 packs per carton

Color: ●



Boot Cover

This Boot Cover offers excellent contamination prevention between footwear and work environment.

Code 2-802

Packing: 100pcs per pack, 10 packs per carton

Color: ●



Plastic Shoe Cover

This Plastic Shoe Cover offers excellent contamination prevention between footwear and work environment.

Code 4-400

Packing: 100pcs per pack, 10 packs per carton

Color: 

Shoe Cover Refill For Mac Shoe Machine

This Plastic Shoe Cover Refill Box for Shoe Cover Machine (MAC-SHOE) offers excellent an contamination prevention between footwear and work environment.

Code MAC-SF

Packing: 100pcs per box, 42 boxes per carton

Color: 

Shoe Cover Machine

This Shoe Cover Machine is automatic machine, easily loaded in refill box of plastic shoe cover, which makes it easy to maintain in higher volume applications. Operated by using electricity, the plastic shoe cover can be easily loaded at this machine. It's ideal for use in a variety of environments such as a laboratory or cleanroom, a manufacturing or food processing plant, a medical or healthcare facility, a film or set location, or anywhere you need a shoe cover dispenser.

Code MAC-SHOE

Packing: 1 unit per box

Dimension: 560x280x380mm


Color: 
Weight: 8.45kg



**TOP
PICKS!**



CROSS PROTECTION[®]



STERILIZATION



Sterilization Pouches

This Pouches are used for packing the medical devices before sterilization and maintain the sterility of the products after sterilization.

CODE	SIZES	PACKING
060130	60mm X 130mm	200 pcs X 30 boxes
070260	70mm X 260mm	200 pcs X 40 boxes
090260	90mm X 260mm	200 pcs X 20 boxes
135280	135mm X 280mm	200 pcs X 20 boxes
135285	135mm X 285mm	200 pcs X 10 boxes
190360	190mm X 360mm	200 pcs X 8 boxes
300490	300mm X 490mm	200 pcs X 5 boxes

Features:

- With triple side for secure closure and double indicator (inside and outside of the pouch).
- CE certified.
- EO and steam (duo indicator).



Sterilization Gussetted Reel

This Gussetted Reel is used for packing the medical devices before sterilization and maintain the sterility of the products after sterilization.

CODE	SIZES	PACKING
G75100	75mm X 38mm X 100m	8 rolls per carton
G10100	100mm X 50mm X 100m	6 rolls per carton
G15100	150mm X 55mm X 100m	4 rolls per carton
G20100	200mm X 65mm X 100m	2 rolls per carton
G25100	250mm X 70mm X 100m	2 rolls per carton
G30100	300mm X 80mm X 100m	2 rolls per carton

Features:

- CE certified.
- EO, Steam, and Form (tri Indicator).



Sterilization Flat Reel

This Flat Reel is used for packing the medical devices before sterilization and maintain the sterility of the products after sterilization.

CODE	SIZES	PACKING
F55200	50mm X 200m	16 rolls per carton
F75200	75mm X 200m	10 rolls per carton
F10200	100mm X 200m	6 rolls per carton
F15200	150mm X 200m	4 rolls per carton
F20200	200mm X 200m	4 rolls per carton
F25200	250mm X 200m	2 rolls per carton
F30200	300mm X 200m	2 rolls per carton
F35200	350mm X 200m	2 rolls per carton
F40200	400mm X 200m	2 rolls per carton

Features:

- CE certified.
- EO, Steam, and Form (tri Indicator).



Sterilization Tyvek Reel

Tyvek Reel is a rather unique material in the stable packaging materials classification. It is known for its outstanding strength, durability and tear resistance. Tyvek Reel is also used for Gas Plasma sterilization at hospitals.

CODE	SIZES	PACKING
P07570	75mm X 70m	8 rolls per carton
P10070	100mm X 70m	6 rolls per carton
P15070	150mm X 70m	4 rolls per carton
P20070	200mm X 70m	4 rolls per carton
P25070	250mm X 70m	2 rolls per carton
P30070	300mm X 70m	2 rolls per carton

Features:

- CE certified.
- Gas plasma Indicator.



Autoclave Tape

Autoclave tapes are manufactured for closing packages of items requiring steam sterilization in hospitals and other similar uses.

CODE	SIZES	PACKING
STR12	12mm X 50m	76 rolls per carton
STR19	19mm X 50m	48 rolls per carton
STR25	25mm X 50m	36 rolls per carton

Green Crepe Paper

Use for packing the medical devices before sterilization and maintain the sterility of the products after sterilization. Suitable for Steam, EO, Plasma

CODE	SIZES	PACKING
CP6060	60cm X 60cm	500 pcs per carton
CP7575	75cm X 75cm	250 pcs per carton
CP9090	90cm X 90cm	250 pcs per carton
CP1010	100cm X 100cm	250 pcs per carton
CP1212	120cm X 120cm	125 pcs per carton
CP1214	120cm X 140cm	125 pcs per carton
CP1216	120cm X 160cm	125 pcs per carton
CP1218	120cm X 180cm	125 pcs per carton

Features:

- Waterproof
- No chips
- Strong bacteria resistance

Sterilization Wrap

CROSS PROTECTION SMMMS Sterilization Wrap is softer and stronger than ever, yet still delivers the barrier protection features to meet your high expectation. Made of medical grade, anti-static material, used to wrap the medical instruments and parts, initial preparation step to go sterilization process.

CODE	SIZES	PACKING
SM5050	50cm X 50cm	500 pcs per carton
SM6060	60cm X 60cm	500 pcs per carton
SM7575	75cm X 75cm	250 pcs per carton
SM9090	90cm X 90cm	250 pcs per carton
SM1010	100cm X 100cm	250 pcs per carton
SM1212	120cm X 120cm	125 pcs per carton
SM1214	120cm X 140cm	125 pcs per carton
SM1216	120cm X 160cm	125 pcs per carton
SM1218	120cm X 180cm	125 pcs per carton

Features:

- Durable sterilization wrap
- Superior strength
- Microbial protection

Material:

- 100% SMMMS

Weight:

- 50GSM



Crepe Bandage

This Crepe Bandage is a woven elasticated bandage and ideal for providing light support to strains and sprains. It is available in a range of width sizes and roll sizes.

5cm x 4.5m

Code 3-200

Packing: 12 rolls per bag, 20 bags per carton

7.5cm x 4.5m

Code 3-201

Packing: 12 rolls per bag, 20 bags per carton

10cm x 4.5m

Code 3-202

Packing: 12 rolls per bag, 20 bags per carton

15cm x 4.5m

Code 3-203

Packing: 12 rolls per bag, 20 bags per carton



Surgical Tape Paper

This Surgical Tape is a type of pressure sensitive adhesive tape used in first aid to hold a bandage or other dressing onto a wound. Our tapes is hypoallergenic adhesive and may remove easily without damaging the skin.

12.5mm x 5m

Code NW12

Packing: 1 pc per card, 360 cards per carton

25mm x 5m

Code NW25

Packing: 1 pc per card, 360 cards per carton

50mm x 5m

Code NW50

Packing: 1 pc per card, 180 cards per carton



Bandage Scissors

Bandage Scissors have large easy grip handles and long curved blades. These all purpose medical scissors are made of stainless steel and can cut thick bandages with ease.

Code 1-124B

Packing: 1 pc per bag, 100 bags per carton

Features:

- Cuts thick bandages
- Easy grip handles
- All purpose scissors
- Great for use at home, travel, or work



Non Woven Sponge Roll (100m)

This Non Woven Sponge roll is highly absorbent and less linting than standard gauze. It is ideal for wound dressing, scrubbing and general purpose cleansing.

Code 3-488

Packing: 1 roll per pack, 12 packs per carton

Dental Roll

This Dental Roll is ideal for used during dental procedure to reduce the saliva flow, keep the mouth open and retract the cheeks.

8x38mm

Code 3-6021

Packing: 500 pcs per bag, 20 bags per carton

10x38mm

Code 3-6022

Packing: 500 pcs per bag, 20 bags per carton

Cotton Wool Roll

400g

Code COT-04

Cotton Wool Roll 400G is mixed with 60% of cotton ideal for first aid, general medical and hygiene.

Packing: 50 rolls per carton

500g

Code 3-601

Cotton Wool Roll 500G is made of 100% cotton and ideal for first aid, general medical and hygiene.

Packing: 12 rolls per carton

Gauze Roll (100yrd)

This Absorbent Gauze Roll made of 100% cotton. Ideal for absorbing fluid as well as dressing and protecting wound. 19x15 mesh, closing weave structure provide excellence absorbency.

Standard - 90cm

Code 3-400

Packing: 1 roll per pack, 12 packs per carton

X-ray Detachable - 90cm

Code 3-401

Packing: 1 roll per pack, 12 packs per carton



Cotton Wool Ball (150gsm)

This 100% pure Cotton Wool Balls are soft, gentle and absorbent. It is perfect for baby and family care, as well as your everyday beauty regime. Each ball is uniformly rolled in consistent shape and size.

Code 3-600

Packing: 250 pcs per bag, 10 bags per carton



Combined Dressing Roll (18cm)

This Dressing Roll is made from 100% pure, non-chlorine bleached cotton fiber. It is ideal for padding or protection for major wounds.

Code 3-489

Packing: 1 roll per pack, 12 packs per carton



Precut Non Woven Sponge 4ply Non Sterile

This Non Sterile Precut Non Woven Sponge is 4ply fold and come with 3 sizes for use in different application.

5cm x 5cm

Code 3-844

Packing: 200 pcs per pack, 50 packs per carton

7.5cm x 7.5cm

Code 3-845

Packing: 200 pcs per pack, 20 packs per carton

10cm x 10cm

Code 3-846

Packing: 200 pcs per pack, 20 packs per carton



Precut Gauze Sponge 8ply Non Sterile

This Non Sterile Precut Gauze Sponge is 8ply fold and come with 3 sizes for use in different application.

5cm x 5cm

Code 3-805

Packing: 200 pcs per pack, 25 packs per carton

7.5cm x 7.5cm

Code 3-875

Packing: 200 pcs per pack, 25 packs per carton

10cm x 10cm

Code 3-810

Packing: 200 pcs per pack, 25 packs per carton



**TOP
PICKS!**

CONSUMER

Uninterrupted sleep &
reduced health risks for
wounds and falls



CAREGIVERS

Improved morale &
productivity

FACILITIES

Reduced product usage,
laundry & odors



CROSS PROTECTION®

COTTON SOFT

Coated Bedsheet &
Pillow Case

**Bacteria Free . Water Proof . Fitted Corner
Breathable . Cloth-Like Back Sheet . Latex Free**

**TOP
PICKS!**

with
**Super Absorbent
Polymer (SAP)**



- ✓ Neutralizes Urine
- ✓ Leak Proff System
- ✓ Bacteria Control
- ✓ Latex-Free



34.0 oz. capacity



CROSS PROTECTION®

Disposable **UNDERPAD**



Single Use Only



Disposable



Non - Sterile

Disposable Underpad, 100% breathable backsheet allows for heat and moisture vapor to escape. Ideal for use with low-air-loss bed systems.

Usage:

- Ideal for moderate to light absorbency protection.
- For routine use as a procedural tray cover or in similar surface protection situation.
- Soft breathable top sheet.
- Absorbent fluff
- Polypropylene backing.

Mini Compressed Towel (White)

This Mini Compressed Towel is great for travel and can use immediate after sop up water.

Code 2-576

Packing: 10 pcs per pack, 100 packs per carton

Size: 30cm x 20cm

Compressed Towel (Assorted Color)

This Compressed Towel is great for travel and can use immediate after sop up water.

Code 2-588

Packing: 10 pcs per bag, 100 bags per carton

Size: 30cm x 60cm

Basic Dressing Set

This non-invasive Dressing Set is ideal for general wound dressing procedure.

Code BAS-1

Packing: 1 set per pouch, 100 pouches per carton

Contain:

ITEMS	QUANTITY
Tray with 3 compartments	1 unit (hard model)
Dressing Forceps	3 unit (grey, green, blue)
Gauze Sponges	4 unit (7.5x7.5cm)
Patient Drape	1 unit (45x45cm)
Sterile Field	1 unit (45x45cm)
Cotton Balls	6 unit (white)
Paper Hand Towel	1 unit (white)
Limpet Bag	1 unit (yellow)

AAMI Disposable Pillow Case

Our AAMI Disposable Pillow Case come with blue color.

Code 2-504

Packing: 10 pcs per bag, 20 bags per carton

Size: 74cm x 48cm

AAMI Disposable Bedsheet

Our AAMI Disposable Bedsheet is fitted corner design and come with blue color.

Code 2-501

Packing: 1 pc per bag, 100 bags per carton

Size: 141cm x 220cm



Cardiology Stethoscope

A stainless steel combination chestpiece for maximum convenience and internal spring binaurals for high acoustic sensitivity. Dual non-chill rings provide patient comfort and latex-free tubing keeps acoustics clear, even in the noisiest settings.

Code 1-125

Packing: 1 pc per box

Features:

- Superior acoustics
- Tunable diaphragm
- Handcrafted solid stainless steel / zinc alloy chestpiece
- Two tubes in one design
- Soft sealing eartips



SAP Absorbent Underpad

Disposable underpads are found in just about any health clinic, hospital or hospice home. The absorbent pads are necessary to keep urine, feces, or any bit of body fluid from damaging skin, soiling bedsheets, staining furniture or flooring. It is not only for human incontinence, but for pets used as pee pads.

17 x 27inch

Code M1727

Packing: 10 pcs per bag, 20 bags per carton

24 x 36inch

Code M2436

Packing: 10 pcs per bag, 12 bags per carton

30 x 30inch

Code M3030

Packing: 10 pcs per bag, 10 bags per carton



Cotton Soft Coated Pillow Case

We are enlisted as a leading manufacturer, supplier and exporter of Pillow Case - Cotton Soft. As per the latest market trends, we have developed our range under the supervisions and able-guidance of the industry heads. These covers are known for their colorfastness, durability and fine finishing.

Code 3-351

Packing: 100pcs per bag

Features:

- Bacteria Free
- Breathable
- Cloth-Like back sheet
- Water Proof
- Latex free

Material: cotton soft

Type: Knitted

Color: ☐



Cotton Soft Coated Bedsheet

We are leading manufacturer and supplier of Bed Sheets - Cotton Soft. We provide the best quality of Hospital bed sheets which are manufactured using high-grade fabrics. These are supplied to various hospitals, nursing homes, and other healthcare centers. Our products are highly durable, comfortable and easy to clean.

Code 3-531

Packing: 100pcs per bag

Features:

- Bacteria Free
- Breathable
- Cloth-Like back sheet
- Fitted Corner
- Water Proof
- Latex Free

Material: Cotton soft

Type: Knitted

Color: ☐



Procedure Set - Cardiovascular Drape

Single use set consists of full sets of components used in cardiovascular procedure. It is intended for use in medical procedure to diagnose and treat some heart conditions.

This set consist of Thoracic Drape, Lower Vascular Drape, Table Cover, Mayo Stand Cover Reinforced, and other types of sets and compartment that are organize in complete for the cardiovascular procedure.

Code PROCEDURE-CARD

Packing: 1 set per pouch

Procedure Set - Caesarean

Caesarean Set provide a sterile field during caesarean procedure. Impervious to strike through, with fluid controls features, reducing the risk of contamination and bacteria transfer.

This set consist of Table Cover, Caesarean Drape, Baby Blanket, Santara Gown & Towel, Electrosurgical Pencil & Tip Cleaner, Umbilical Cord Lamp, Suction Catheter & Tubing, Syringe, Urine Bag, Foley Catheter, and other types of sets and compartment that are organize in complete for the caesarean procedure.

Code CAESAREAN-CP

Packing: 1 set per pouch

Procedure Set - Ophthalmic

For use cataract removal, corneal graft, retinal detachment or basic eye procedure.

This set consist of Surgical Gown, Eye Drape, Gauze, Syringe, Cotton Applicator, and other types of sets and compartment that are organize in complete for the ophthalmic procedure.

Code OPH-01

Packing: 1 set per pouch



Dental Mirror

This Dental Mirror is one of the most common instruments used in the dental clinics, a part of diagnostic instruments in Dentistry.

Code 5-102

Packing: 100 pcs per box, 10 boxes per carton

Color: ○

Dental Probe

This Dental Probe is primary used to measure pocket depths around a tooth in order to establish the state of health of the periodontium.

Code 5-101S

Packing: 100 pcs per box, 10 boxes per carton

Color: ●

Microtip Surgical Aspirator (Sterile)

To facilitate the connection between the surgical aspirator & the dental suction unit.

Code 100-2367

Packing: 25 pcs per box, 10 boxes per carton

Denture Box

This Denture Box is Ideal for keeping your denture clean and safe.

Code A777

Packing: 100 pcs per carton

Color: ●●●●●

Dappen Dish

Our Dappen Dish is made of high quality plastic which is impervious to all dental materials. The dishes are designed for quick and easy mixing of materials such as prophylactic paste, fluoride products, composite materials, acrylics, cements and amalgams.

Code 5-200

Packing: 50 pcs per bag, 20 bags per carton

Color: ●●●●●○



TOP
PICKS!



CROSS PROTECTION®

SALIVA
EJECTOR



Saliva Ejector

Our Saliva Ejector 145mm is made of certified medical and food grade PVC. The selected raw materials, patented design featuring round edges, smooth surfaces and perfectly placed suction ports that help to enhance patient comfort while protecting delicate mucosal tissue.

Code 5-501

Packing: 100 pcs per bag, 10 bags per carton

Size: 145mm

Features:

- Maximum Safety
- Smooth Edge
- Bend Easily
- Smooth Suction, No Block

TIP COLOR	TUBE COLOR
Assorted color	Clear
Blue	Blue
Blue	Clear
White	White
White	Clear
Green	Green
Pink	Pink
Yellow	Yellow

Dental Tubing Sleeve

Tubing sleeve helps to prevent contamination and reduces the need for dangerous, harsh chemical sanitation through its disposability. Disposable anti-reflective transparent polyethylene tubing sleeve with open ends. Easy to use.

Code 550-182

Packing: 1 roll per box, 10 boxes per carton

Size: 2" diameter x 100" length per roll

Dental Bib

Made of 2 ply paper and 1 ply PE. Extra heavy weight 18-18gsm origin paper provides excellent absorbency.

Code 5-800

Packing: 50 pcs per bag, 10 bags per carton

Size: 45cm x 33cm

Color: 

Barrier Film Dispenser

Clear

Durable plastic film dispenser to organize your barrier films. Available in clear acrylic.

Code 5-388D

Packing: 1 pc per box

White

Elegant and strong barrier film dispenser that can fit regular size wrapping rolls.

Code 5-388H

Packing: 1 pc per box



Paper Hand Towel

Our Paper Hand Towel is designed for clean up spills and moisture. It is suitable for general usage.

C Fold

Code M1753

Packing: 180 pcs per bag, 20 bags per carton

Color: ☐

Size: 23x30cm

Feature:

- High Absorbency
- Durable
- Extra Thickness
- High White

M Fold

Code M1753-M

Packing: 250 pcs per bag, 12 bags per carton

Color: ☐

Size: 23x21cm

Feature:

- High Absorbency
- Durable
- Extra Thickness
- High White

Hand Roll Towel

Our Hand Roll Towel is suitable for hand drying.

Code M1843

Packing: 8 rolls per carton

Color: ☐

Size: 200 meters

Feature:

- High Absorbency
- Durable
- Extra Thickness
- High White

Protective Goggles Clear Antifog

Protective Goggles - Clear Antifog feature large ventilated vinyl frames that can be worn over regular glasses and have easily adjusted straps for comfort and security.

Code 5-606

Packing: 1pcs per pack, 20 packs per carton



Paper Hand Towel Dispenser

This M Fold Paper Hand Towel Dispenser dispenses folded paper towel for hand drying and wiping up spills. Paper can easily be taken from the unit without having to touch it, promoting good hygienic habits.

Code M1743-DM

Packing: 1 unit per box

Hand Roll Towel Dispenser (AUTOCUT)

This Autocut Hand Roll Towel Dispenser is recommended to use at office, hotel, hospital, retail shop, F&B, factory, shopping mall and etc.

Code M1843-D

Packing: 1 unit per box

Waste Bag

Small

Code	456515B	Color	Blue
Code	456515L	Color	Black

Packing: 30 pcs per bag, 60 bags per carton
Size: 450 x 650 x 0.015mm

Medium

Code	658502B	Color	Blue
------	---------	-------	------

Packing: 20 pcs per bag, 50 bags per carton
Size: 650 x 850 x 0.02mm

Large

Code	801002B	Color	Blue
------	---------	-------	------

Packing: 10 pcs per bag, 30 bags per carton
Size: 800 x 1000 x 0.025mm

X-Large

Code	9011003B	Color	Blue
Code	9011003L	Color	Black

Packing: 5 pcs per bag, 60 bags per carton
Size: 900 x 1100 x 0.03mm



First Aid Bag Small Compact

Code FIRST-01

Packing: 100 sets per carton

Size: 13cm x 20cm

Contains:

ITEMS	QUANTITY
Bandage 5cmx4.5m	1
Face Mask (Ear Loop)	1
Antiseptic Cleansing Wipes	4
Triangular Bandage 96x96x135cm	1
Surgical Tape	1
Alcohol Pad	6
First Aid Scissors	1
Safety Pins	3
Adhesive Strips	10
Tourniquet	1
Wound Pad (Sterile)	2
Cotton Ball (Sterile) 0.5g/ball	10
Forceps	1
Zip Lock Bag	1
First Aid Bag	1



First Aid Bag Large

Code FIRST-02

Packing: 50 sets per carton

Size: 16cm x 45cm

Contains:

ITEMS	QUANTITY
Disposable Gloves	1
Cool Fever	1
Zip Lock Bag	1
Bandage 75cmx4.5cm	2
Face Mask (Ear Loop)	1
Sting Relief Pad	4
Triangular Bandage 96x96x135cm	1
Surgical Tape	1
Tourniquet	1
First Aid Scissors	1
Safety Pins	3
Non Woven Dressing (6x8cm)	2
Quick Saver (Mouth-mouth life saving device)	1
Wounds Pad (Sterile)	4
Antiseptic Cleansing Wipes	2
Forceps	1
Alcohol Pad	6
Cotton Ball (Sterile) 0.5g/ball	10
Thermometer	1
Adhesive Strips	10
First Aid Bag	1



Shop **ONLINE** with us

- ✓ Sign up today
- ✓ Convenience
- ✓ Save up to 15% and more



Sign up /
Login



Browse
Products



Add to
Cart



Payment



Shipping
Methods



Checkout

www.crossprotection.com.my

Tel: +603-3392 0000

Fax: +603-3392 0001



CROSS PROTECTION (M) SDN BHD
No.11, Jalan Perindustrian 5, Batu 5 1/2,
Off Jalan Haji Abdul Manan, 41050 Klang,
Selangor, Malaysia.



operation@crossprotection.com /
marketing@crossprotection.com /
sales@crossprotection.com



+603 33920000 / +6016 2620436

www.crossprotection.com.my